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
Dear customer, thank you for purchasing our product. Please read the operating instructions carefully before putting the appliance into operation and keep these instructions including the warranty, the receipt and, if possible, the box with the internal packing.

General Safety Instructions


- The appliance is designed exclusively for private use and for the envisaged purpose.
- This appliance is not fit for commercial use.
- Keep it away from sources of heat, direct sunlight, humidity (never dip it into any liquid) and sharp edges. Do not use the appliance with wet hands.
- If you should leave the workplace, always switch the machine off.
- Do not try to repair the appliance on your own. Always contact an authorized technician.
- Pay careful attention to the following “Special Safety Instructions for this unit“.

Special safety instructions for this unit


- Operate the appliance only in accordance with the specifications on the rating plate.
- Do not put the appliance on the wet places.
- Clean the scale with a slightly damp cloth. Do not immerse the scale in the water or use chemical / abrasive agents.
- Always use the scael on a hard, flat surface. DO NOT use on carpet
- DO NOT use the scale on wet and slippery floor. DO NOT step on the edge of the scale. The scale platform is slippery when wet. Keep it dry!.
- Always weigh yourself with the same scales places on the same surface. Uneven floors may affect the reading.
- Weigh yourself without clothing and footwear, before meals and always at the same time of day.
- Pacemaker users / wearers, or other medical implants with electronic components, are not recommended to use this scale.
- The body fat & hydration% are for reference only. Please always consult your doctor for advice of medication or diet in order to achieve ideal body fat level.
- To make sure accuracy of testing, please DO NOT use the scale in environment with strong electro-magnetic field.
- DO NOT use the scale on wet and slippery floor. DO NOT step on the edge of the scale. The scale platform is slippery when wet. Keep it dry!
- DO NOT use this scale on a place with intense shake. Keep stand still on scale when testing.
- Do not drop the device, or jump on the device when using. Avoid strong impact and sudden shake as it is a precision measuring device. The construction is delicate and sensitive to all unwanted external impact.
- This appliance is not intended for use by persons (including children) with reduced physical. Sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
- Children should be supervised to ensure that they do not play with the appliance.

Description of the controls

P1. Unit of scale
P2. Display
P3. Glass surface
P4. Electrodes
P5. ON / SET switch ON/SET = ■
P6. START switch START = ▶
P7. UP and DOWN arrow press switch ▼ ▲

LCD segmentum / icons illustrations:

Specifications:

Max. Weighing Capacity: 150 kg, 330 lb
Min. Weighing at: 2.0 kg, 4,4 lb
Weight Measuring Division: 100 g, 0,22 lb
Body Fat measure range: 4.0-45.0%
Hydration measure range: 37.8-66.0%
Hydration differential resolution in %: 0.1%

Parameter	Default setting	Range
Memory	1	1 - 8
Gender	Male	Male/ Female
Height	165 cm	100 - 250 cm
Age	25	10 - 100

Instructions for Use

Introduction Of Body Fat Testing
1. About Body Fat Percentage
Body fat percentage indicates body fat proportion in body, which is a criteria to judge one is overweight or not.
2. Testing Principle
Impedance of fat is higher than any other body tissues such as muscle, when electrical signal of a certain frequency goes through body (BIA principle) . Based on this principle, body fat percentage can be tested by sending an electrical signal through body.

The technology used in these scales, Bio-electrical Impedance Analysis (BIA) is more technologically advanced than the order Body Mass Index (BMI) method where fat % is calculated using only height and weight .BIA technology uses 5 parameters to get a more accurate resolution for the analysis (weight and biological resistance which are automatically measured by the scales together with the three personal parameters entered manually - sex, height and age).

Important:

During body fat testing, please remove shoes and socks, etc., and clear your feet, in order to get the most accurate result.
To ensure accuracy, please test at the same time of the day
Test result may be inaccurate when the users are:
- afte r intense sport activity, overfeeding, or excessive dehydration.
- children under 10, or adults over 80;
- body building atheletes, or other occupational atheletes;
- pregnant women; Dialysis or edema patient.

Normal weighing scale operation

1) Turn on the unit by slight tap on the lower right corner of the platform of the scale. LCD will wake up prompt to General Weighing Mode, showing „0.0 kg“ (or „0.0 lb“)
2) Step on the scale gently (Do not jump on the scale). As soon as you step on the scale, LCD will immediately display incrementing digits till it stabilises followed by flashes and then locked digits, showing your weight value. After a few seconds, LCD will automatically shut down.

Fat monitor scale operation

1/ Press button (■) Key to wake up the unit. All LCD segmentum light up briefly and prompt to setting mode starting from Memory selection. By pressing button (■) key again will take you through cycle of setting all parameters. See following illustration:

a) Press Up / Down Arrow (▲▼) button to advance / reverse to the memory station which you would like to choose or has been designated to you. Press (■) button to confirm setting. LCD will prompt to next parameter setting (Gender) at the same time.

b) Press Up / Down Arrow (▲▼) button to advance / reverse to choose the gender of your type. Press (■) button to confirm setting. LCD will prompt to next parameter setting (Height) at the same time.

c) Press Up / Down Arrow (▲▼) button to advance / reverse to the adjust to your height level. Press (■) button to confirm setting, LCD will prompt to next parameter (Age) at the same time.

d) Press Up / Down Arrow (▲▼) button to advance/reverse to adjust to your age. Press (■) button to confirm setting. LCD will prompt to Memory setting. All settings complete!!! At this point, you can continue to press (■) button a few consecutive times to re-cycle through the setting mode to review and assure all the parameters are correctly set.

e) Press (▶) button to exit setting mode and enter normal Fat monitor operating mode. LCD will prompt to ‚0.0 kg‘ (or ‚0.0‘ lb). You are ready to weigh yourself, and let the scale compute your Fat %, Hydration %.

2/ Step on the scale gently (ensure to remove your socks & shoes , wipe clean your feet & the weighing platform to get the best conductive effect). While stepping on, align your feet to the metal conductive pads, two for each foot. Ensure your feet are standing apart not touching each other (to avoid short circuit)

If you step on the scale too late after the LCD auto-off, simply press the (■) button again, wait till LCD prompt to show your current settings. Then press (▶) button again to re-enter normal fat monitor scale operation mode.

3/ As soon as you step on the scale, LCD will immediately display incrementing digits till it stabilises, locked & flashes, indicating your current weight.

4/ Within a few seconds, the scale will automatically compute your Fat % and Hydration %. LCD will prompt to show your Fat % (icon „BF“) followed by hydration % (icon „BW“). Subsequently, LCD will prompt back and forth alternately between value of Fat % and hydration % for 5 times, and then automatically shut down.

By the lower edge of the LCD, you will see a key line segments of multiple vertical bars. This key line serves to give an indication of your Fat % level, telling where your fat level standing is comparing to the statistics norms. Whether you are considered as underweight, normal, or overweight. Your can take this indication as a reference data of your condition.

NOTE:
Whenever pressing (■) button, turning on the unit for fresh operation, you can always go through the same steps from Illustration (a) to (e) to verify or adjust all personal data and previous measured records, before measuring for new records

Error and other LCD message

LCD may show occasionally error message‘ Err ,. Simply turn off the unit for a brief moment, or remove and re-install battery as a total re-set. Then turn on again to resume operation.

1. All segments lit up unanimously: Indicating all segments performing correctly.
2. „ - - - - - „ = Strain gauge sensor internal zero point setting too low. Can not operate properly.
3. „ Err „= Overload, or strain gauge sensor internal zero point too high .
4. „ LO „ = Main unit battery power running low.
5. „ Err-h „ = The computed Fat % rate too high above range of upper limit.
6. „ Err-L „ = The computed Fat % rate too low below range of lower limit.

Cleaning and Maintenance

1) Use a dry, soft fabric to clean the device. Or a lightly damp soft fabric with mild detergent to wipe clean the device. Do not use strong chemical solution such as alcohol, benzene for cleaning.
2) It could be slippery when the unit platform surface is wet, or semi-damp. Ensure the weighing platform surface is dry before using for safety reason.

Troubleshooting

The device is not working
- Check the rights position of battery and the polarity.
- When warning indicator shows „Lo“ replace the battery.


Technical specifications

- Electronic body fat and water scale
- Capacity 150 kg, 330 lb
- Graduation 100 g, 0,22 lb
- Measurement units kg/lb
- Large display 76 x 38 mm
- 8 mm tempered glass surface
- Using 5 parameters for measuring of body fat – height, weight, age, sex and biological resistant
- Memory for 8 person
- Foot tap switch
- Auto off, auto zero
- Overload indicator
- Low battery indicator
- Batteries operated 9 V (not included)
- Size: 35 x 4,2 x 32 cm
- Weight: 2,44 kg

We reserve the right to change technical specifications.

WARNING: Do not use this product near water, in wet areas to avoid fire or injury of electric current. Always turn off the product when you don't use it or before a revision. There aren't any parts in this appliance which are reparable by consumer. Always appeal to a qualified authorized service. The product is under a dangerous tention.

Liquidation of old electrical and electronic devices

This symbol appearing on the product, on the product accessories or on the product packing means that the product must not be disposed as household waste. Please transport the product to the respective collection point where the electric and electronic device will be recycled. Collection points for the used electric and electronic devices exist in the European Union and also in other European countries. By proper liquidation of the product you can prevent possible negative impacts on environment and human health, which might otherwise occur as a consequence of improper manipulation with the product. Recycling of materials contributes to protection of natural resources. Therefore please do not throw the old electric and electronic devices in the household wastes. For more detailed information please contact your local authority, any service company liquidating household waste or the store where you have bought the product.

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